

For Research Use Only

WDR19 Monoclonal antibody, PBS Only



Catalog Number: 67330-1-PBS

Basic Information

Catalog Number:

67330-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG4599

GenBank Accession Number:

BC032578

GeneID (NCBI):

57728

UNIPROT ID:

Q8NEZ3

Full Name:

WD repeat domain 19

Calculated MW:

1342 aa, 152 kDa

Observed MW:

145 kDa

Purification Method:

Protein A purification

CloneNo.:

2D1H4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67330-1-Ig (WDR19 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67330-1-PBS in a different storage buffer formulation.