

For Research Use Only

PSMD9 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67338-3-PBS



Basic Information

Catalog Number: 67338-3-PBS	GenBank Accession Number: BC004213	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 5715	CloneNo.: 1D5C2
Isotype: IgG1	UNIPROT ID: O00233	
Immunogen Catalog Number: AG25654	Full Name: proteasome (prosome, macropain) 26S subunit, non-ATPase, 9	
	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51695-1, PSMD9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67338-2-PBS. Detection antibody: 67338-3-PBS.



Cytometric bead array standard curve of MP51695-1, PSMD9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67338-2-PBS. Detection antibody: 67338-3-PBS. Standard: Ag25654. Range: 0.781-100 ng/mL.