

For Research Use Only

N-PAC Monoclonal antibody

Catalog Number: 67340-1-Ig



Basic Information

Catalog Number:

67340-1-Ig

Size:

1700 μ g/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG7133

GenBank Accession Number:

BC064940

GeneID (NCBI):

84656

UNIPROT ID:

Q49A26

Full Name:

cytokine-like nuclear factor n-pac

Calculated MW:

61 kDa

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

1A5E1

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Mouse, Rat

Positive Controls:

WB : HeLa cells, HEK-293 cells, HepG2 cells, Jurkat cells, HSC-T6 cells, NIH/3T3 cells, 4T1 cells

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67340-1-Ig (N-PAC antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.