

For Research Use Only

N-PAC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67340-2-PBS



Basic Information

Catalog Number:

67340-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7133

GenBank Accession Number:

BC064940

GeneID (NCBI):

84656

UNIPROT ID:

Q49A26

Full Name:

cytokine-like nuclear factor n-pac

Calculated MW:

61 kDa

Purification Method:

Protein G purification

CloneNo.:

1G8G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51700-1, N-PAC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67340-2-PBS. Detection antibody: 67340-3-PBS.



Cytometric bead array standard curve of MP51700-1, N-PAC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67340-2-PBS. Detection antibody: 67340-3-PBS. Standard: Ag7133. Range: 1.563-200 ng/mL.