

For Research Use Only

N-PAC Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67340-3-PBS



Basic Information

Catalog Number: 67340-3-PBS	GenBank Accession Number: BC064940	Purification Method: Protein A Magarose purification
Source: Mouse	GeneID (NCBI): 84656	CloneNo.: 2B1D2
Isotype: IgG2a	UNIPROT ID: Q49A26	
Immunogen Catalog Number: AG7133	Full Name: cytokine-like nuclear factor n-pac	
	Calculated MW: 61 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51700-1, N-PAC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67340-2-PBS. Detection antibody: 67340-3-PBS.



Cytometric bead array standard curve of MP51700-1, N-PAC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67340-2-PBS. Detection antibody: 67340-3-PBS. Standard: Ag7133. Range: 1.563-200 ng/mL.