

For Research Use Only

ATG7 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67341-3-PBS



Basic Information

Catalog Number:

67341-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29172

GenBank Accession Number:

BC000091

GeneID (NCBI):

10533

UNIPROT ID:

O95352

Full Name:

ATG7 autophagy related 7 homolog
(S. cerevisiae)

Calculated MW:

78 kDa

Purification Method:

Protein G purification

CloneNo.:

3E6F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51585-1, ATG7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67341-3-PBS. Detection antibody: 67341-4-PBS.



Cytometric bead array standard curve of MP51585-1, ATG7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67341-3-PBS. Detection antibody: 67341-4-PBS. Standard: Ag29172. Range: 1.563-200 ng/mL.