

For Research Use Only

TRIM2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67342-2-PBS



Basic Information

Catalog Number:

67342-2-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG14637

GenBank Accession Number:

BC011052

GeneID (NCBI):

23321

UNIPROT ID:

Q9C040

Full Name:

tripartite motif-containing 2

Calculated MW:

744 aa, 82 kDa

Purification Method:

Protein A purification

CloneNo.:

2C9C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51520-1, TRIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67342-1-PBS. Detection antibody: 67342-2-PBS.



Cytometric bead array standard curve of MP51520-1, TRIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67342-1-PBS. Detection antibody: 67342-2-PBS. Standard:Ag14637. Range: 12.5-400 ng/mL.