

For Research Use Only

MUS81 Monoclonal antibody, PBS Only



Catalog Number: 67351-1-PBS

Featured Product

Basic Information

Catalog Number:

67351-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG29089

GenBank Accession Number:

BC009999

GeneID (NCBI):

80198

UNIPROT ID:

Q96NY9

Full Name:

MUS81 endonuclease homolog (S. cerevisiae)

Calculated MW:

61 kDa

Observed MW:

52-72 kDa

Purification Method:

Protein A purification

CloneNo.:

2D9B7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67351-1-Ig (MUS81 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67351-1-PBS in a different storage buffer formulation.



WB result of MUS81 antibody (67351-1-Ig; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-MUS81 transfected MCF-7 cells. This data was developed using the same antibody clone with 67351-1-PBS in a different storage buffer formulation.