

For Research Use Only

# ARL4C Monoclonal antibody

Catalog Number: 67356-1-Ig



## Basic Information

Catalog Number:

67356-1-Ig

Size:

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29496

GenBank Accession Number:

BC001051

GeneID (NCBI):

10123

UNIPROT ID:

P56559

Full Name:

ADP-ribosylation factor-like 4C

Calculated MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

1B8B5

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

