

For Research Use Only

RALBP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67359-3-PBS



Basic Information

Catalog Number: 67359-3-PBS	GenBank Accession Number: BC013126	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 10928	CloneNo.: 2A7C9
Isotype: IgG1	UNIPROT ID: Q15311	
Immunogen Catalog Number: AG28618	Full Name: raIA binding protein 1	
	Calculated MW: 655 aa, 76 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51591-1, RALBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67359-2-PBS. Detection antibody: 67359-3-PBS.

Cytometric bead array standard curve of MP51591-1, RALBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67359-2-PBS. Detection antibody: 67359-3-PBS. Standard: Ag28618. Range: 3.125-400 ng/mL.