

For Research Use Only

GAS6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67367-2-PBS



Basic Information

Catalog Number:

67367-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5158

GenBank Accession Number:

BC038984

GeneID (NCBI):

2621

UNIPROT ID:

Q14393

Full Name:

growth arrest-specific 6

Calculated MW:

678 aa, 75 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D5E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50908-1, GAS6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67367-2-PBS. Detection antibody: 67367-3-PBS. Standard:Ag5158. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50908-2, GAS6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67367-2-PBS. Detection antibody: 67367-4-PBS. Standard:Ag5158. Range: 0.098-3.125 ng/mL.