

For Research Use Only

# BRD4 Monoclonal antibody

Catalog Number: 67374-1-Ig **2 Publications**



## Basic Information

Catalog Number:

67374-1-Ig

Size:

1000  $\mu$ g/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27545

GenBank Accession Number:

BC035266

GeneID (NCBI):

23476

UNIPROT ID:

O60885

Full Name:

bromodomain containing 4

Calculated MW:

1362 aa, 152 kDa

Observed MW:

150 kDa

Purification Method:

Protein A purification

CloneNo.:

1F11A6

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

Human

Cited Species:

human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yanwen Lan	32596881	J Cell Mol Med	WB,IF
Zhixin Qiu	34615656	Cancer Discov	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

