

For Research Use Only

# BRD4 Monoclonal antibody, PBS Only



Catalog Number: 67374-2-PBS

## Basic Information

Catalog Number:

67374-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27545

GenBank Accession Number:

BC035266

GeneID (NCBI):

23476

UNIPROT ID:

O60885

Full Name:

bromodomain containing 4

Calculated MW:

1362 aa, 152 kDa

Observed MW:

200 kDa

Purification Method:

Protein G purification

CloneNo.:

2F7E2

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67374-2-Ig (BRD4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67374-2-PBS in a different storage buffer formulation.



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 67374-2-Ig (BRD4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67374-2-PBS in a different storage buffer formulation.