

For Research Use Only

HNRNPH1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67375-2-PBS



Basic Information

Catalog Number:

67375-2-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6603

GenBank Accession Number:

BC001348

GeneID (NCBI):

3187

UNIPROT ID:

P31943

Full Name:

heterogeneous nuclear
ribonucleoprotein H1 (H)

Calculated MW:

49 kDa

Purification Method:

Protein A purification

CloneNo.:

3G11A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51734-1, HNRNP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67375-1-PBS. Detection antibody: 67375-2-PBS.



Cytometric bead array standard curve of MP51734-1, HNRNP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67375-1-PBS. Detection antibody: 67375-2-PBS. Standard: Ag6603. Range: 3.125-400 ng/mL.