

For Research Use Only

ZO-2 Monoclonal antibody, PBS Only

Catalog Number: 67387-1-PBS

Featured Product



Basic Information

Catalog Number:

67387-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG13577

GenBank Accession Number:

BC027592

GeneID (NCBI):

9414

UNIPROT ID:

Q9UDY2

Full Name:

tight junction protein 2 (zona occludens 2)

Calculated MW:

134 kDa

Observed MW:

150 kDa

Purification Method:

Protein A purification

CloneNo.:

1D2E5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67387-1-Ig (ZO-2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67387-1-PBS in a different storage buffer formulation.



WB result of ZO-2 antibody (67387-1-Ig; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-ZO-2 transfected HeLa cells. This data was developed using the same antibody clone with 67387-1-PBS in a different storage buffer formulation.