

For Research Use Only

RGS14 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67394-3-PBS



Basic Information

Catalog Number:

67394-3-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9477

GenBank Accession Number:

BC014094

GeneID (NCBI):

10636

UNIPROT ID:

O43566

Full Name:

regulator of G-protein signaling 14

Calculated MW:

566 aa, 61 kDa

Purification Method:

Protein A purification

CloneNo.:

2F7F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51623-2, RGS14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67394-2-PBS. Detection antibody: 67394-3-PBS.



Cytometric bead array standard curve of MP51623-2, RGS14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67394-2-PBS. Detection antibody: 67394-3-PBS. Standard:Ag9477. Range: 3.125-400 ng/mL