

For Research Use Only

CYP21A2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67421-2-PBS



Basic Information

Catalog Number:

67421-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26292

GenBank Accession Number:

NM_000500

GeneID (NCBI):

1589

UNIPROT ID:

P08686

Full Name:

cytochrome P450, family 21,
subfamily A, polypeptide 2

Calculated MW:

56 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D5A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50770-1, CYP21A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67421-2-PBS. Detection antibody: 67421-3-PBS. Standard:Ag26292. Range: 3.125-100 ng/mL.