

For Research Use Only

# EIF4G2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67428-2-PBS



## Basic Information

Catalog Number:

67428-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12080

GenBank Accession Number:

BC043149

GeneID (NCBI):

1982

UNIPROT ID:

P78344

Full Name:

eukaryotic translation initiation factor 4 gamma, 2

Calculated MW:

907 aa, 102 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A3E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP51146-1, EIF4G2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67428-2-PBS. Detection antibody: 67428-3-PBS. Standard:Ag12080. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP51146-3, EIF4G2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67428-2-PBS. Detection antibody: 67428-5-PBS. Standard:Ag12080. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP51146-4, EIF4G2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67428-2-PBS. Detection antibody: 67428-6-PBS. Standard:Ag12080. Range: 3.125-100 ng/mL