

For Research Use Only

EIF4G2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67428-3-PBS



Basic Information

Catalog Number:	67428-3-PBS	GenBank Accession Number:	BC043149	Purification Method:	Protein G purification
Size:	1 mg/ml	GeneID (NCBI):	1982	CloneNo.:	2F2B4
Source:	Mouse	UNIPROT ID:	P78344		
Isotype:	IgG1	Full Name:	eukaryotic translation initiation factor 4 gamma, 2		
Immunogen Catalog Number:	AG12080	Calculated MW:	907 aa, 102 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51146-1, EIF4G2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67428-2-PBS. Detection antibody: 67428-3-PBS. Standard:Ag12080. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP51146-2, EIF4G2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67428-4-PBS. Detection antibody: 67428-3-PBS. Standard:Ag12080. Range: 0.781-100 ng/mL.