

For Research Use Only

EIF4G2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67428-5-PBS



Basic Information

Catalog Number: 67428-5-PBS	GenBank Accession Number: BC043149	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1982	CloneNo.: 2A1H6
Source: Mouse	UNIPROT ID: P78344	
Isotype: IgG2a	Full Name: eukaryotic translation initiation factor 4 gamma, 2	
Immunogen Catalog Number: AG12080	Calculated MW: 907 aa, 102 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51146-3, EIF4G2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67428-2-PBS. Detection antibody: 67428-5-PBS. Standard:Ag12080. Range: 0.781-100 ng/mL.