

For Research Use Only

# EIF4G2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67428-6-PBS



## Basic Information

Catalog Number:	67428-6-PBS	GenBank Accession Number:	BC043149	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	1982	CloneNo.:	2H5B1
Source:	Mouse	UNIPROT ID:	P78344		
Isotype:	IgG2b	Full Name:	eukaryotic translation initiation factor 4 gamma, 2		
Immunogen Catalog Number:	AG12080	Calculated MW:	907 aa, 102 kDa		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP51146-4, EIF4G2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67428-2-PBS. Detection antibody: 67428-6-PBS. Standard:Ag12080. Range: 3.125-100 ng/mL.