

For Research Use Only

# NUDT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67443-2-PBS



## Basic Information

Catalog Number:

67443-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9802

GenBank Accession Number:

BC014618

GeneID (NCBI):

4521

UNIPROT ID:

P36639

Full Name:

nudix (nucleoside diphosphate linked moiety X)-type motif 1

Calculated MW:

179aa, 20 kDa; 197aa, 23 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D6G1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50783-1, NUDT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67443-2-PBS. Detection antibody: 67443-3-PBS. Standard:Ag9802. Range: 0.195-50 ng/mL