

For Research Use Only

# HNRNPA2B1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67445-4-PBS



## Basic Information

Catalog Number:

67445-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6699

GenBank Accession Number:

BC000506

GeneID (NCBI):

3181

UNIPROT ID:

P22626

Full Name:

heterogeneous nuclear  
ribonucleoprotein A2/B1

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

3G5E11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51320-2, HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67445-2-PBS. Detection antibody: 67445-4-PBS.



Cytometric bead array standard curve of MP51320-2, HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67445-2-PBS. Detection antibody: 67445-4-PBS. Standard:Ag6699. Range: 0.195-100 ng/mL.