

For Research Use Only

# MYC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67447-2-PBS



## Basic Information

Catalog Number:

67447-2-PBS

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17519

GenBank Accession Number:

BC000141

GeneID (NCBI):

4609

UNIPROT ID:

P01106

Full Name:

v-myc myelocytomatosis viral  
oncogene homolog (avian)

Calculated MW:

49 kDa

Purification Method:

Protein G purification

CloneNo.:

2E12E12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50347-1, MYC/c-MYC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67447-2-PBS. Detection antibody: 67447-3-PBS. Standard: Ag17519. Range: 0.781-50 ng/mL