

For Research Use Only

MYC Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67447-3-PBS



Basic Information

Catalog Number:

67447-3-PBS

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17519

GenBank Accession Number:

BC000141

GeneID (NCBI):

4609

UNIPROT ID:

P01106

Full Name:

v-myc myelocytomatosis viral oncogene homolog (avian)

Calculated MW:

49 kDa

Purification Method:

Protein G purification

CloneNo.:

3D2A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50347-1, MYC/c-MYC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67447-2-PBS. Detection antibody: 67447-3-PBS. Standard:Ag17519. Range: 0.781-50 ng/mL.