

For Research Use Only

HSP90AB1 Monoclonal antibody

Catalog Number: 67450-1-Ig

Featured Product

1 Publications



Basic Information

Catalog Number:

67450-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28379

GenBank Accession Number:

BC014485

GeneID (NCBI):

3326

UNIPROT ID:

P08238

Full Name:

heat shock protein 90kDa alpha
(cytosolic), class B member 1

Calculated MW:

723 aa, 83 kDa

Observed MW:

83-90 kDa

Purification Method:

Protein A purification

CloneNo.:

2D3D1

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells, HEK-293 cells, HepG2 cells, Jurkat cells, HSC-T6 cells, NIH/3T3 cells, Neuro-2a cells

Background Information

HSP90AB1, also known as HSP90 beta, is a member of HSP90 family which includes cytosolic HSP90 alpha (HSP90AA1) and HSP90 beta, endoplasmic reticulum related GRP74, mitochondrial TRAP1. HSP90 proteins play a versatile role in cell regulation, forming complexes with a large number of cellular kinases, transcription factors, and other molecules. Upregulation of HSP90AB1 has been found in various cancers.

Notable Publications

Author	Pubmed ID	Journal	Application
Bo Wang	38008382	Virol Sin	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67450-1-Ig (HSP90AB1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



WB result of HSP90AB1 antibody (67450-1-Ig; 1:20000; incubated at room temperature for 1.5 hours) with sh-Control and sh-HSP90AB1 transfected HEK-293 cells.