

For Research Use Only

WNT3 Monoclonal antibody

Catalog Number: 67452-1-Ig **4 Publications**



Basic Information

Catalog Number:

67452-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG27957

GenBank Accession Number:

BC114219

GeneID (NCBI):

7473

UNIPROT ID:

P56703

Full Name:

wingless-type MMTV integration site family, member 3

Calculated MW:

355 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

2B4D11

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

Human, mouse, rat, pig

Cited Species:

human, mouse

Positive Controls:

WB : pig brain tissue, mouse brain tissue, rat brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yunhe An	32536825	Cancer Cell Int	WB
Meng-Xuan Gui	33420580	Chin J Integr Med	IHC, WB
Xin-Zhan Zhang	38994173	World J Gastrointest Oncol	

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67452-1-Ig (WNT3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.