

For Research Use Only

AXIN2 Monoclonal antibody, PBS Only

Catalog Number: 67458-1-PBS



Basic Information

Catalog Number:

67458-1-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29291

GenBank Accession Number:

BC006295

GeneID (NCBI):

8313

UNIPROT ID:

Q9Y2T1

Full Name:

axin 2

Calculated MW:

843 aa, 94 kDa

Observed MW:

100 kDa

Purification Method:

Protein G purification

CloneNo.:

1A11A4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

Axis inhibition protein2 (AXIN2), also known as Cductin or Axil, is a multidomain scaffold protein that negatively regulate Wnt signaling. AXIN2 can directly interact with beta-catenin and GSK3B. AXIN2 has a number of phosphorylation sites, and also can undergo poly(ADP-ribose)lation by tankyrase TNKS and TNKS2. Poly(ADP-ribose)lated AXIN2 then get ubiquitinated by RNF146 and leads to its degradation and subsequent activation of Wnt signaling. AXIN2 is localized in the cytoplasm and its mutation is involved in colorectal cancer.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67458-1-Ig (AXIN2 antibody) at dilution of 1:4800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67458-1-PBS in a different storage buffer formulation.