

For Research Use Only

# Arfaptin-1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67461-3-PBS



## Basic Information

Catalog Number:

67461-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG12097

GenBank Accession Number:

BC103759

GeneID (NCBI):

27236

UNIPROT ID:

P53367

Full Name:

ADP-ribosylation factor interacting protein 1

Calculated MW:

373 aa, 42 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E3H3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51228-1, Arfaptin-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67461-2-PBS. Detection antibody: 67461-3-PBS. Standard:Ag12097. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51228-2, Arfaptin-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67461-3-PBS. Detection antibody: 67461-1-PBS. Standard:Ag12097. Range: 3.125-100 ng/mL