

For Research Use Only

SLC25A12 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67467-4-PBS



Basic Information

Catalog Number:

67467-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13576

GenBank Accession Number:

BC016932

GeneID (NCBI):

8604

UNIPROT ID:

O75746

Full Name:

solute carrier family 25
(mitochondrial carrier, Aralar),
member 12

Calculated MW:

75 kDa, 63 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E3B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50650-2, SLC25A12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67467-4-PBS. Detection antibody: 67467-1-PBS. Standard:Ag13576. Range: 1.563-100 ng/mL