

For Research Use Only

# SF3B3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67469-3-PBS



## Basic Information

Catalog Number:

67469-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6315

GenBank Accession Number:

BC003146

GeneID (NCBI):

23450

UNIPROT ID:

Q15393

Full Name:

splicing factor 3b, subunit 3, 130kDa

Calculated MW:

136 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3H9E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51555-1, SF3B3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67469-2-PBS. Detection antibody: 67469-3-PBS.



Cytometric bead array standard curve of MP51555-1, SF3B3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67469-2-PBS. Detection antibody: 67469-3-PBS. Standard:Ag6315. Range: 3.125-400 ng/mL.