

For Research Use Only

# FURIN Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67481-3-PBS



## Basic Information

Catalog Number:

67481-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13313

GenBank Accession Number:

BC012181

GeneID (NCBI):

5045

UNIPROT ID:

P09958

Full Name:

furin (paired basic amino acid cleaving enzyme)

Calculated MW:

87 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G5D4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51279-1, FURIN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67481-2-PBS. Detection antibody: 67481-3-PBS. Standard:Ag13313. Range: 1.563-100 ng/mL.