

For Research Use Only

S100P Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67485-2-PBS



Basic Information

| | | | | | |
|---------------------------|-------------|---------------------------|--------------------------------|----------------------|------------------------|
| Catalog Number: | 67485-2-PBS | GenBank Accession Number: | BC006819 | Purification Method: | Protein G purification |
| Size: | 1 mg/ml | GeneID (NCBI): | 6286 | CloneNo.: | 2E7G2 |
| Source: | Mouse | UNIPROT ID: | P25815 | | |
| Isotype: | IgG1 | Full Name: | S100 calcium binding protein P | | |
| Immunogen Catalog Number: | AG19115 | Calculated MW: | 95 aa, 11 kDa | | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51168-1, S100P Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67485-2-PBS. Detection antibody: 67485-3-PBS. Standard:Ag19115. Range: 0.781-100 ng/mL.