

For Research Use Only

# RHOBTB3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67502-3-PBS



## Basic Information

Catalog Number:

67502-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5782

GenBank Accession Number:

BC041337

GeneID (NCBI):

22836

UNIPROT ID:

O94955

Full Name:

Rho-related BTB domain containing 3

Calculated MW:

611 aa, 69 kDa

Purification Method:

Protein G purification

CloneNo.:

1G4F9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51247-1, RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67502-2-PBS. Detection antibody: 67502-3-PBS. Standard:Ag5782. Range: 0.391-100 ng/mL