

For Research Use Only

RHOBTB3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67502-5-PBS



Basic Information

Catalog Number:

67502-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5782

GenBank Accession Number:

BC041337

GeneID (NCBI):

22836

UNIPROT ID:

O94955

Full Name:

Rho-related BTB domain containing 3

Calculated MW:

611 aa, 69 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B12A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51247-3, RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67502-5-PBS. Detection antibody: 67502-6-PBS. Standard:Ag5782. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP51247-4, RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67502-5-PBS. Detection antibody: 67502-7-PBS. Standard:Ag5782. Range: 1.563-100 ng/mL