

For Research Use Only

CDK5R1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67503-3-PBS



Basic Information

Catalog Number: 67503-3-PBS	GenBank Accession Number: BC020580	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 8851	CloneNo.: 2G3A9
Isotype: IgG1	UNIPROT ID: Q15078	
Immunogen Catalog Number: AG13152	Full Name: cyclin-dependent kinase 5, regulatory subunit 1 (p35)	
	Calculated MW: 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51665-1, CDK5R1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67503-2-PBS. Detection antibody: 67503-3-PBS.



Cytometric bead array standard curve of MP51665-1, CDK5R1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67503-2-PBS. Detection antibody: 67503-3-PBS. Standard: Ag13152. Range: 1.563-200 ng/mL.