

For Research Use Only

SIRT6 Monoclonal antibody, PBS Only

Catalog Number: 67510-1-PBS

Featured Product



Basic Information

Catalog Number:

67510-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24103

GenBank Accession Number:

BC028220

GeneID (NCBI):

51548

UNIPROT ID:

Q8N6T7

Full Name:

sirtuin (silent mating type
information regulation 2 homolog) 6
(S. cerevisiae)

Calculated MW:

39 kDa

Observed MW:

40 kDa

Purification Method:

Protein G purification

CloneNo.:

2H4G8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

SIRT6, also named as SIR2L6, belongs to the sirtuin family. It is known to deacetylate H3K9, H3K56 and CtIP in vivo. SIRT6 regulates the expression of a large number of stress-responsive and metabolism related genes, promotes genomic stability and stimulates base excision and double strand break DNA repair(PMID: 20729089). By maintaining both the integrity and the expression of the mammalian genome, SIRT6 thus serves several roles that parallel Sir2 function(PMID: 20729089). In addition, Sirt6-deficient animals develop normally for several weeks after birth, but succumb to an acute multi-organ degenerative syndrome that is uniformly lethal by one month of age(PMID: 19135889).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67510-1-Ig (SIRT6 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67510-1-PBS in a different storage buffer formulation.



WB result of SIRT6 antibody (67510-1-Ig; 1:4000; incubated at room temperature for 1.5 hours) with sh-Control and sh-SIRT6 transfected Jurkat cells. This data was developed using the same antibody clone with 67510-1-PBS in a different storage buffer formulation.