

For Research Use Only

BSEP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67512-2-PBS



Basic Information

Catalog Number: 67512-2-PBS	GenBank Accession Number: NM_003742	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 8647	CloneNo.: 2C7C1
Isotype: IgG1	UNIPROT ID: O95342	
Immunogen Catalog Number: AG29135	Full Name: ATP-binding cassette, sub-family B (MDR/TAP), member 11	
	Calculated MW: 146 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

