

For Research Use Only

BSEP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67512-3-PBS



Basic Information

Catalog Number:

67512-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29135

GenBank Accession Number:

NM_003742

GeneID (NCBI):

8647

UNIPROT ID:

O95342

Full Name:

ATP-binding cassette, sub-family B (MDR/TAP), member 11

Calculated MW:

146 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B3E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

