

For Research Use Only

FES Monoclonal antibody, PBS Only

Catalog Number: 67524-1-PBS



Basic Information

Catalog Number:

67524-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28108

GenBank Accession Number:

BC035357

GeneID (NCBI):

2242

UNIPROT ID:

P07332

Full Name:

feline sarcoma oncogene

Calculated MW:

822 aa, 94 kDa

Observed MW:

93 kDa, 79 kDa, 87 kDa, 85 kDa

Purification Method:

Protein A purification

CloneNo.:

1E1B11

Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67524-1-Ig (FES antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67524-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using 67524-1-Ig (FES antibody), at dilution of 1:300 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67524-1-PBS in a different storage buffer formulation.