

For Research Use Only

GSDMB Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67527-2-PBS



Basic Information

Catalog Number:

67527-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28332

GenBank Accession Number:

BC025682

GeneID (NCBI):

55876

UNIPROT ID:

Q8TAX9

Full Name:

gasdermin B

Calculated MW:

411 aa, 47 kDa

Purification Method:

Protein G purification

CloneNo.:

1B12G9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array sample test of MP51638-1, GSDMB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67527-2-PBS. Detection antibody: 67527-3-PBS.



Cytometric bead array standard curve of MP51638-1, GSDMB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67527-2-PBS. Detection antibody: 67527-3-PBS. Standard:Ag28332. Range: 3.125-400 ng/mL.



Bi-layer interferometry (BLI) kinetic assays of 67527-2-Ig against Human GSDMB were performed. The affinity constant is 25.6 nM.