

For Research Use Only

Tim23 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67535-4-PBS



Basic Information

Catalog Number:

67535-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27069

GenBank Accession Number:

BC062707

GeneID (NCBI):

100287932

UNIPROT ID:

O14925

Full Name:

translocase of inner mitochondrial membrane 23

Calculated MW:

22 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3F4D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50554-2, Tim23 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67535-3-PBS. Detection antibody: 67535-4-PBS. Standard:Ag27069. Range: 0.098-100 ng/mL.