

For Research Use Only

KLRB1/CD161 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67537-3-PBS



Basic Information

Catalog Number:

67537-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG27117

GenBank Accession Number:

BC114516

GeneID (NCBI):

3820

ENSEMBL Gene ID:

ENSG00000111796

UNIPROT ID:

Q12918

Full Name:

killer cell lectin-like receptor
subfamily B, member 1

Calculated MW:

225 aa, 25 kDa

Purification Method:

Euglobulin precipitation

CloneNo.:

3F10A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50517-1, KLRB1,CD161 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67537-2-PBS. Detection antibody: 67537-3-PBS. Standard:Ag27117. Range: 0.781-25 ng/mL.