

For Research Use Only

RBM8A,Y14 Monoclonal antibody, PBS Only (Capture)

Catalog Number:67541-2-PBS



Basic Information

Catalog Number: 67541-2-PBS	GenBank Accession Number: BC017088	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 9939	CloneNo.: 2B11D1
Isotype: IgG1	UNIPROT ID: Q9Y5S9	
Immunogen Catalog Number: AG6800	Full Name: RNA binding motif protein 8A	
	Calculated MW: 20 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51528-1, RBM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67541-2-PBS. Detection antibody: 67541-3-PBS.



Cytometric bead array standard curve of MP51528-1, RBM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67541-2-PBS. Detection antibody: 67541-3-PBS. Standard:Ag6800. Range: 1.563-200 ng/mL