

For Research Use Only

RBM8A,Y14 Monoclonal antibody, PBS Only (Detector)

Catalog Number:67541-3-PBS



Basic Information

Catalog Number:

67541-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6800

GenBank Accession Number:

BC017088

GeneID (NCBI):

9939

UNIPROT ID:

Q9Y5S9

Full Name:

RNA binding motif protein 8A

Calculated MW:

20 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A5B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51528-1, RBM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67541-2-PBS. Detection antibody: 67541-3-PBS.



Cytometric bead array standard curve of MP51528-1, RBM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67541-2-PBS. Detection antibody: 67541-3-PBS. Standard: Ag6800. Range: 1.563-200 ng/mL.