

For Research Use Only

CRX Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67555-2-PBS



Basic Information

Catalog Number: 67555-2-PBS	GenBank Accession Number: BC016664	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 1406	CloneNo.: 1C8G2
Source: Mouse	UNIPROT ID: O43186	
Isotype: IgG1	Full Name: cone-rod homeobox	
Immunogen Catalog Number: AG30086	Calculated MW: 299 aa, 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51402-1, CRX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67555-2-PBS. Detection antibody: 67555-3-PBS.



Cytometric bead array standard curve of MP51402-1, CRX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67555-2-PBS. Detection antibody: 67555-3-PBS. Standard: Ag30086. Range: 0.098-100 ng/mL.