

For Research Use Only

BTN3A3 Monoclonal antibody

Catalog Number: 67560-1-Ig



Basic Information

Catalog Number:

67560-1-Ig

Size:

1700 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12755

GenBank Accession Number:

BC015815

GeneID (NCBI):

10384

UNIPROT ID:

O00478

Full Name:

butyrophilin, subfamily 3, member A3

Calculated MW:

584 aa, 65 kDa

Observed MW:

65 kDa

Purification Method:

Protein G purification

CloneNo.:

1A3B5

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, HepG2 cells, Jurkat cells, K-562 cells, HL-60 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67560-1-Ig (BTN3A3 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.