

For Research Use Only

# TAF6L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67569-3-PBS



## Basic Information

Catalog Number:

67569-3-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG8202

GenBank Accession Number:

BC008785

GeneID (NCBI):

10629

UNIPROT ID:

Q9Y6J9

Full Name:

TAF6-like RNA polymerase II,  
p300/CBP-associated factor (PCAF)-  
associated factor, 65kDa

Calculated MW:

622 aa, 68 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F10E9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51658-1, TAF6L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67569-2-PBS. Detection antibody: 67569-3-PBS.



Cytometric bead array standard curve of MP51658-1, TAF6L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67569-2-PBS. Detection antibody: 67569-3-PBS. Standard: Ag8180. Range: 1.563-200 ng/mL.