

For Research Use Only

TBC1D25 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67574-3-PBS



Basic Information

Catalog Number: 67574-3-PBS	GenBank Accession Number: BC101817	Purification Method: Protein A Magarose purification
Source: Mouse	GeneID (NCBI): 4943	CloneNo.: 3F10H1
Isotype: IgG2b	UNIPROT ID: Q3MII6	
Immunogen Catalog Number: AG29470	Full Name: TBC1 domain family, member 25	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51606-1, TBC1D25 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67574-2-PBS. Detection antibody: 67574-3-PBS.



Cytometric bead array standard curve of MP51606-1, TBC1D25 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67574-2-PBS. Detection antibody: 67574-3-PBS. Standard: Ag29470. Range: 0.781-100 ng/mL.