

For Research Use Only

CDK1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67575-2-PBS



Basic Information

Catalog Number:

67575-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29018

GenBank Accession Number:

BC014563

GeneID (NCBI):

983

UNIPROT ID:

P06493

Full Name:

cell division cycle 2, G1 to S and G2 to M

Calculated MW:

25 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F9C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51605-1, CDK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67575-2-PBS. Detection antibody: 67575-3-PBS.



Cytometric bead array standard curve of MP51605-1, CDK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67575-2-PBS. Detection antibody: 67575-3-PBS. Standard: Ag29018. Range: 0.781-100 ng/mL.